

Anti-Cyp3A4 Rabbit Polyclonal Antibody For Wb And Elisa - P08684

Item Number: CM0015500



Introduction

Rabbit polyclonal antibody targeting human CYP3A4, validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for studying xenobiotic metabolism, steroid hormone biosynthesis, and drug metabolism research.

[Learn More](#)

Parameter	Value
Product Name	CYP3A4 Rabbit pAb
Remarks/Alias	HLP; CP33; CP34; CYP3A; NF-25; VDDR3; CYP3A3; P450C3; CYP11A3; CYP11A4; P450PCN1; CYP3A4
Species	Human
GeneID (human)	1576
GeneID	1576
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 244-503 of human CYP3A4 (NP_059488.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08684
Protein Weight	57kDa
Shipping	Ice bag