


Anti-Cyp3A5 Polyclonal Antibody For Wb, Elisa - P20815

Item Number: CM0027043



Introduction

Rabbit polyclonal antibody targeting human CYP3A5, validated for WB and ELISA. Recognizes human, mouse, and rat CYP3A5. Ideal for studying drug metabolism and steroid hormone pathways. - P20815

[Learn More](#)

Parameter	Value
Product Name	CYP3A5 Rabbit pAb
Remarks/Alias	CP35; PCN3; CYP3A5; P450PCN3; CYP3A5
Species	Human
Gene ID	1577 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 233-502 of human CYP3A5 (NP_000768.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20815
Protein Weight	57kDa
Shipping	Ice bag