

Anti-Cyp46A1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9Y6A2

Item Number: CM0022494



Introduction

CYP46A1 Rabbit Polyclonal Antibody validated for WB, IF-P, ELISA; reacts with human, mouse, rat. Targets cholesterol 24-hydroxylase, a brain cholesterol homeostasis enzyme. Immunogen: synthetic peptide within aa 400-500. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	CYP46A1 Rabbit pAb
Remarks/Alias	CP46; CYP46; CYP46A1
Species	Human
GeneID	10858
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-500 of human CYP46A1 (NP_006659.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y6A2
Protein Weight	57kDa
Shipping	Ice bag