


Anti-Cyp46A1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y6A2

Item Number: CM0027776



Introduction

Rabbit polyclonal antibody against human CYP46A1, a cholesterol 24-hydroxylase. Validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for neuroscience research on cholesterol metabolism in the brain.

[Learn More](#)

Parameter	Value
Product Name	CYP46A1 Rabbit pAb
Remarks/Alias	CP46, CYP46, CYP46A1
Species	Human
GeneID (Human)	10858
GeneID	10858
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-280 of human CYP46A1 (NP_006659.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6A2
Protein Weight	57kDa
Shipping	Ice bag