

Anti-Cyp4F11 Polyclonal Antibody For If/Icc, Elisa - Q9Hbi6

Item Number: CM0023774

Introduction

Rabbit polyclonal antibody against human CYP4F11, validated for IF/ICC and ELISA. Targets a cytochrome P450 enzyme involved in fatty acid metabolism, vitamin K inactivation, and drug demethylation. Suitable for liver and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	CYP4F11 Rabbit pAb
Remarks/Alias	CYP4F11; CYP4F11
Species	Human
GeneID (Human)	57834
GeneID	57834
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-130 of human CYP4F11 (NP_001122404.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9HBI6
Protein Weight	60kDa
Shipping	Ice bag