

Anti-Cyp4F12 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Hcs2

Item Number: CM0025630

Introduction

Anti-CYP4F12 rabbit polyclonal antibody for WB and ELISA. Cross-reacts with human, mouse, and rat. Targets cytochrome P450 4F12 involved in fatty acid and drug metabolism.

[Learn More](#)

Parameter	Value
Product Name	CYP4F12 Rabbit pAb
Remarks/Alias	CYP4F12; F22329_1; CYP4F12
Species	Human
GeneID (human)	66002
GeneID	66002
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 265-524 of human CYP4F12 (NP_076433.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HCS2
Protein Weight	60kDa
Shipping	Ice bag