

Anti-Cyp51A1 Polyclonal Antibody For Wb, Elisa - Q16850

Item Number: CM0029428



Introduction

Rabbit polyclonal antibody targeting CYP51A1 (lanosterol 14- α demethylase), validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for cholesterol biosynthesis and sterol metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	CYP51A1 Rabbit pAb
Remarks/Alias	LDM; CP51; CYP51; CYPL1; P450L1; P450-14DM; CYP51A1
Species	Human
GeneID (Human)	1595
GeneID	1595
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 210-509 of CYP51A1 (NP_000777.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16850
Protein Weight	57kDa
Shipping	Ice bag