


Anti-Cyp51A1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16850

Item Number: CM0013522



Introduction

CYP51A1 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Detects lanosterol 14-alpha demethylase, a cholesterol biosynthesis enzyme. Ideal for metabolism research.

[Learn More](#)

Parameter	Value
Product Name	CYP51A1 Rabbit pAb
Remarks/Alias	LDM; CP51; CYP51; CYPL1; P450L1; P450-14DM; CYP51A1
Species	Human
Gene ID (Human)	1595
Gene ID	1595
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 245-420 of CYP51A1 (NP_000777.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16850
Protein Weight	57kDa
Shipping	Ice bag