


Anti-Cyp51A1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q16850

Item Number: CM0007592



Introduction

Rabbit monoclonal antibody against human CYP51A1 (lanosterol 14-alpha demethylase), raised against recombinant protein, for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Critical for cholesterol biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	CYP51A1 Rabbit mAb
Remarks/Alias	LDM; CP51; CYP51; CYPL1; P450L1; P450-14DM; CYP51A1
Species	Human
Gene ID (Human)	1595
Gene ID	1595
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16850
Protein Weight	57kDa
Shipping	Ice bag