

Anti-Cyp7A1 Rabbit Polyclonal Antibody For If-P And Elisa - P22680

Item Number: CM0026564



Introduction

CYP7A1 Rabbit Polyclonal Antibody validated for IF-P and ELISA in human samples. Targets cholesterol 7- α -hydroxylase, a key enzyme in bile acid biosynthesis and cholesterol homeostasis, essential for metabolism research.

[Learn More](#)

Parameter	Value
Product Name	CYP7A1 Rabbit pAb
Remarks/Alias	CP7A; CYP7; CYPVII
Species	Human
GeneID (Human)	1581
GeneID	1581
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CYP7A1 (NP_000771.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Human
SWISS	P22680
Protein Weight	58kDa
Shipping	Ice bag