

Anti-Cyp7B1 Polyclonal Antibody For Wb, If/Icc, Elisa - 075881

Item Number: CM0018796

Introduction

Rabbit polyclonal antibody targeting human CYP7B1, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for oxysterol metabolism, neurosteroid synthesis, and B cell migration research. - 075881

[Learn More](#)

Parameter	Value
Product Name	CYP7B1 Rabbit pAb
Remarks/Alias	CP7B; CBAS3; SPG5A; CYP7B1
Species	Human
Gene ID (Human)	9420
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CYP7B1 (NP_004811.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75881
Protein Weight	58kDa
Shipping	Ice bag