

Anti-Cypor Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P16435

Item Number: CM0011957



Introduction

Rabbit anti-CYPOR monoclonal antibody for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets NADPH--cytochrome P450 reductase, essential for electron transfer from NADP.

[Learn More](#)

Parameter	Value
Product Name	CYPOR Rabbit mAb
Remarks/Alias	CPR; CYPOR; P450R
Species	Human
GeneID (Human)	5447
GeneID	5447
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human CYPOR (P16435).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16435
Protein Weight	77 kDa
Shipping	Ice bag