

Anti-Cypor Rabbit Polyclonal Antibody For Wb And Elisa - P16435

Item Number: CM0021681



Introduction

Rabbit polyclonal antibody against human CYPOR (NADPH--cytochrome P450 reductase), validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studying electron transfer in microsomal cytochrome P450 systems.

[Learn More](#)

Parameter	Value
Product Name	CYPOR Rabbit pAb
Remarks/Alias	CPR; CYPOR; P450R
Species	Human
GeneID (Human)	5447
GeneID	5447
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human CYPOR (NP_000932.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16435
Protein Weight	77kDa
Shipping	Ice bag