

Anti-Cyr61 Polyclonal Antibody For Wb, If/Icc, Elisa - 000622

Item Number: CM0017733



Introduction

Rabbit polyclonal antibody against human CYR61 (CCN1/GIG1/IGFBP10) for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. 42kDa secreted protein involved in cell adhesion and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	CYR61 Rabbit pAb
Remarks/Alias	GIG1; CYR61; IGFBP10
Species	Human
GeneID (Human)	3491
GeneID	3491
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 120-220 of human CYR61 (NP_001545.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00622
Protein Weight	42kDa
Shipping	Ice bag