


Anti-Cysltr1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y271

Item Number: CM0022977



Introduction

Rabbit polyclonal antibody against cysteinyl leukotriene receptor 1 (CYSLTR1), validated for WB and ELISA in human and mouse. Suitable for studying inflammation, asthma, and GPCR signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	CYSLTR1 Rabbit pAb
Remarks/Alias	CYSLT1; CYSLTR; CYSLT1R; HMTMF81; CYSLTR1
Species	Human
GeneID (Human)	10800
GeneID	10800
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human CYSLTR1 (NP_006630.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y271
Protein Weight	39kDa
Shipping	Ice bag