


## Anti-Cystathionase/Cth Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P32929

Item Number: CM0025115



### Introduction

Cystathionase/CTH Rabbit pAb is a polyclonal antibody validated for WB, IF/ICC, IP, ELISA in human, mouse, rat. It targets cystathionine gamma-lyase, a key enzyme in trans-sulfuration and H<sub>2</sub>S production.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | Cystathionase/CTH Rabbit pAb                                                                                          |
| Remarks/Alias    | CGL; CSE; Cystathionase/CTH                                                                                           |
| Species          | Human                                                                                                                 |
| GeneID (Human)   | 1491                                                                                                                  |
| Immunogen        | Recombinant protein containing a sequence corresponding to amino acids 1-405 of human Cystathionase/CTH (NP_001893.2) |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, IF/ICC, IP, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                     |
| SWISS            | P32929                                                                                                                |
| Protein Weight   | 45kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |