

Anti-Cystathionine Beta-Synthase (Cbs) Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P35520

Item Number: CM0002757



Introduction

Rabbit monoclonal antibody targeting human cystathionine beta-synthase (CBS, P35520), validated for WB, IHC-P, IF/ICC, ELISA. Suitable for research on homocysteine metabolism, transsulfuration pathway, and hydrogen sulfide signaling.

[Learn More](#)

Parameter	Value
Product Name	CBS Rabbit mAb
Remarks/Alias	CBSL; HIP4; CBS
Species	Human
GeneID (Human)	875
GeneID	875
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-551 of human CBS (P35520)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P35520
Protein Weight	61kDa
Shipping	Ice bag