

Anti-Cystatin A Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P01040

Item Number: CM0010107



Introduction

Rabbit monoclonal antibody against human Cystatin A (CSTA). Validated for WB, IHC-P, IF/ICC, ELISA, and cross-reacts with mouse and rat. Ideal for studying protease inhibition and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	Cystatin A Rabbit mAb
Remarks/Alias	AREI; PSS4; STF1; STFA; Cystatin A
Species	Human
Gene ID (Human)	1475
Gene ID	1475
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-98 of human Cystatin A (NP_005204.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01040
Protein Weight	11kDa
Shipping	Ice bag