

Anti-Cystatin B Monoclonal Antibody For Wb, Elisa - P04080

Item Number: CM0003967



Introduction

High-quality rabbit monoclonal antibody against human Cystatin B (Stefin-B), validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studying protease inhibition and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	Cystatin B Rabbit mAb
Remarks/Alias	PME; ULD; CST6; EPM1; STFB; CPI-B; EPM1A; Cystatin B
Species	Human
Gene ID (Human)	1476
Gene ID	1476
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-98 of human Cystatin B (NP_000091.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04080
Protein Weight	11kDa
Shipping	Ice bag