


Anti-Cystatin C Polyclonal Antibody For Wb, Elisa - P01034

Item Number: CM0013123



Introduction

High-quality rabbit anti-Cystatin C polyclonal antibody validated for Western blot and ELISA. Reacts with human, mouse, and rat Cystatin C, a secreted cysteine protease inhibitor implicated in amyloidosis and age-related macular degeneration.

[Learn More](#)

Parameter	Specification
Product Name	Cystatin C Rabbit pAb
Remarks/Alias	ARMD11; HEL-S-2; Cystatin C
Species	Human
GeneID (Human)	1471
GeneID	1471
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-146 of human Cystatin C (NP_000090.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01034
Protein Weight	13kDa
Shipping	Ice bag