


Anti-Cystatin C Rabbit Monoclonal Antibody For Wb, Elisa - P01034

Item Number: CM0012586



Introduction

High-quality rabbit monoclonal anti-Cystatin C antibody (Clone A8933, P01034) validated for WB and ELISA. Detects human and mouse cystatin C; ideal for renal and neurological disease research.

[Learn More](#)

Parameter	Value
Product Name	Cystatin C Rabbit mAb
Remarks/Alias	ARMD11; HEL-S-2; Cystatin C
Species	Human
Gene ID (Human)	1471
Gene ID	1471
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-146 of human Cystatin C (P01034)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01034
Protein Weight	13kDa
Shipping	Ice bag