


Anti-Cystatin C Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - P01034

Item Number: CM0004964



Introduction

Recombinant rabbit monoclonal antibody against human Cystatin C, validated for WB, IF, IHC, and ELISA. Cross-reacts with human and mouse. Ideal for studying cysteine protease regulation and amyloidosis.

[Learn More](#)

Parameter	Value
Product Name	Cystatin C Rabbit mAb
Remarks/Alias	ARMD11, HEL-S-2, Cystatin C
Species	Human
Gene ID	1471
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-146 of human Cystatin C (NP_000090.1)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01034
Protein Weight	13kDa
Shipping	Ice bag