

Anti-Cytochrome C Monoclonal Antibody For Wb, Ihc-P, Elisa - P99999

Item Number: CM0010307



Introduction

Human Cytochrome C (P99999) rabbit monoclonal antibody, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial electron transport and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Cytochrome C Rabbit mAb
Remarks/Alias	CYC; HCS; THC4; C
Species	Human
GeneID (Human)	54205
GeneID	54205
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-105 of human Cytochrome C (P99999)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P99999
Protein Weight	12kDa
Shipping	Ice bag