

Anti-Cytochrome C Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P99999

Item Number: CM0013459



Introduction

Rabbit polyclonal antibody targeting human Cytochrome C, validated for WB, IHC-P, IF/ICC, and ELISA applications. Detects human, mouse, and rat Cytochrome C. Ideal for studying apoptosis and mitochondrial electron transport.

[Learn More](#)

Parameter	Value
Product Name	Cytochrome C Rabbit pAb
Remarks/Alias	CYC; HCS; THC4; C
Species	Human
GeneID (Human)	54205
GeneID	54205
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human Cytochrome C (NP_061820.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P99999
Protein Weight	12kDa
Shipping	Ice bag