

## Anti-Cytochrome C Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P99999

Item Number: CM0015360



### Introduction

CamelBio's anti-Cytochrome C rabbit polyclonal antibody is validated for WB, IF/ICC, IHC-P, and ELISA. It detects human, mouse, and rat Cytochrome C, a key apoptosis and mitochondrial respiration marker.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Cytochrome C Rabbit pAb                                                                                                 |
| Remarks/Alias    | CYC; HCS; THC4; Cytochrome C                                                                                            |
| Species          | Human                                                                                                                   |
| Gene ID (Human)  | 54205                                                                                                                   |
| Gene ID          | 54205                                                                                                                   |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human Cytochrome C (NP_061820.1) |
| Source           | Rabbit                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                   |
| Application      | WB, IF/ICC, IHC-P, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                       |
| SWISS            | P99999                                                                                                                  |
| Protein Weight   | 12 kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                 |