

Anti-Cytochrome P450 2E1 (Cyp2E1) Rabbit Mab For Wb, Ihc-P, Elisa - P05181

Item Number: CM0007152



Introduction

Rabbit monoclonal antibody against human CYP2E1, validated in WB, IHC-P, ELISA; cross-reacts with mouse and rat. Suitable for studying fatty acid metabolism and xenobiotic detoxification.

[Learn More](#)

Parameter	Value
Product Name	Cytochrome P450 2E1 (CYP2E1) Rabbit mAb
Remarks/Alias	CPE1; CYP2E; P450-J; P450C2E; Cytochrome P450 2E1 (CYP2E1)
Species	Human
GeneID (Human)	1571
GeneID	1571
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 234-493 of human CYP2E1 (NP_000764.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05181
Protein Weight	57kDa
Shipping	Ice bag