

Anti-Cytohesin-1 (Cyth1) Polyclonal Antibody For Wb, If/Icc, Elisa - Q15438

Item Number: CM0021032



Introduction

Rabbit polyclonal antibody targeting human Cytohesin-1 (CYTH1). Validated for WB, IF/ICC, ELISA, cross-reacts with mouse and rat. ~46 kDa protein. Ideal for ARF activation, membrane trafficking, and epithelial polarization studies.

[Learn More](#)

Parameter	Value
Product Name	CYTH1 Rabbit pAb
Remarks/Alias	B2-1; SEC7; PSCD1; D17S811E; CYTOHESIN-1; CYTH1
Species	Human
GeneID (Human)	9267
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human CYTH1 (NP_004753.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15438
Protein Weight	46kDa
Shipping	Ice bag