


Anti-Cytokeratin 1 (Krt1) Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P04264

Item Number: CM0004514



Introduction

High-quality anti-Cytokeratin 1 (KRT1) rabbit monoclonal antibody validated for WB, IF-P, and ELISA applications. Reacts with human, mouse, and rat. Ideal for studies on keratinocyte differentiation and epidermal structural integrity.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 1 (KRT1) Rabbit mAb
Remarks/Alias	K1; CK1; EHK; AEI2; EHK1; EPPK; KRT1A; NEPPK; Cytokeratin 1 (KRT1)
Species	Human
Gene ID (Human)	3848
Gene ID	3848
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 545-644 of human Cytokeratin 1 (KRT1) (P04264)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04264
Protein Weight	66kDa
Shipping	Ice bag