


Anti-Cytokeratin 10 Monoclonal Antibody For If-P, Elisa - P13645

Item Number: CM0004493



Introduction

Recombinant rabbit monoclonal antibody targeting human Cytokeratin 10 (CK-10, K10). Validated for immunofluorescence (IF-P) and ELISA, cross-reacts with mouse and rat. Supports keratinocyte differentiation and skin barrier research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 10 Rabbit mAb
Remarks/Alias	BIE; EHK; IHL; K10; KPP; BCIE; CK10; Cytokeratin 10
Species	Human
GeneID (Human)	3858
GeneID	3858
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 485-584 of human Cytokeratin 10 (P13645).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13645
Protein Weight	59kDa
Shipping	Ice bag