

## Anti-Cytokeratin 12 (Krt12) Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - Q99456

Item Number: CM0002480

### Introduction

Rabbit monoclonal antibody targeting human Cytokeratin 12 (KRT12), validated for IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Essential for studies of corneal epithelium, intermediate filaments, and Meesmann corneal dystrophy.

[Learn More](#)

| Parameter        | Value                                                                                                                               |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Cytokeratin 12 (KRT12) Rabbit mAb                                                                                                   |
| Remarks/Alias    | K12; MECD1; Cytokeratin 12 (KRT12)                                                                                                  |
| Species          | Human                                                                                                                               |
| Gene ID (Human)  | 3859                                                                                                                                |
| Gene ID          | 3859                                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 250-494 of human Cytokeratin 12 (KRT12) (NP_000214.1) |
| Source           | Rabbit                                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                                               |
| Application      | IF-P, IHC-P, ELISA                                                                                                                  |
| Cross Reactivity | Mouse, Rat                                                                                                                          |
| SWISS            | Q99456                                                                                                                              |
| Protein Weight   | 54kDa                                                                                                                               |
| Shipping         | Ice bag                                                                                                                             |