

Anti-Cytokeratin 12 (Krt12) Rabbit Polyclonal Antibody For Wb And Elisa - Q99456

Item Number: CM0029441

Introduction

High-quality rabbit polyclonal antibody targeting human Cytokeratin 12 (KRT12). Suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for corneal epithelium research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 12 (KRT12) Rabbit pAb
Remarks/Alias	K12; MECD1; Cytokeratin 12 (KRT12)
Species	Human
GeneID (human)	3859
GeneID	3859
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 295-494 of human Cytokeratin 12 (Cytokeratin 12 (KRT12)) (NP_000214.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99456
Protein Weight	54kDa
Shipping	Ice bag