

Anti-Cytokeratin 13 (Krt13) Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P13646

Item Number: CM0002117

Introduction

Rabbit monoclonal antibody against human Cytokeratin 13 (KRT13), validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for mucosal epithelial research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 13 (KRT13) Rabbit mAb
Remarks/Alias	K13; CK13; WSN2; Cytokeratin 13 (KRT13)
Species	Human
Gene ID (Human)	3860
Gene ID	3860
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13646
Protein Weight	50 kDa
Shipping	Ice bag