

Anti-Cytokeratin 13 (Krt13) Rabbit Polyclonal Antibody For Wb, Elisa - P13646

Item Number: CM0028071

Introduction

Rabbit polyclonal antibody targeting human Cytokeratin 13 (KRT13, CK13, K13), validated for Western blot and ELISA. Suitable for human, mouse, rat samples. MW ~50kDa. Research use in mucosal epithelial biology.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 13 (KRT13) Rabbit pAb
Remarks/Alias	K13; CK13; WSN2; Cytokeratin 13 (KRT13)
Species	Human
Gene ID (Human)	3860
Gene ID	3860
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 269-458 of human Cytokeratin 13 (Cytokeratin 13 (KRT13)) (NP_705694.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13646
Protein Weight	50kDa
Shipping	Ice bag