

Anti-Cytokeratin 13 (Krt13) Rabbit Polyclonal Antibody For Wb, Ihc-P, If-P, Elisa - P13646

Item Number: CM0020846

Introduction

Rabbit polyclonal antibody against human Cytokeratin 13 (KRT13), also known as CK13 or WSN2. Suitable for Western blot, IHC, IF, and ELISA applications in epithelial and mucosal research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 13 (KRT13) Rabbit pAb
Remarks/Alias	K13; CK13; WSN2; Cytokeratin 13 (KRT13)
Species	Human
GeneID (Human)	3860
GeneID	3860
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 269-458 of human Cytokeratin 13 (KRT13) (NP_705694.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P13646
Protein Weight	50kDa
Shipping	Ice bag