


Anti-Cytokeratin 14 (Krt14) Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P02533

Item Number: CM0018682



Introduction

Rabbit polyclonal antibody against human Cytokeratin 14 (KRT14/CK14) validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, and rat KRT14. Ideal for research on keratinization, epithelial integrity, and epidermolysis bullosa.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 14 (KRT14) Rabbit pAb
Remarks/Alias	K14; NFJ; CK14; EBS1; EBS3; EBS4; EBS1A; EBS1B; EBS1C; EBS1D; Cytokeratin 14 (KRT14)
Species	Human
Gene ID (human)	3861
Gene ID	3861
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 371-470 of human Cytokeratin 14 (KRT14) (NP_000517.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02533
Protein Weight	52kDa
Shipping	Ice bag