


Anti-Cytokeratin 14 (Krt14) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P02533

Item Number: CM0003344



Introduction

Cytokeratin 14 (KRT14) Rabbit Monoclonal Antibody (P02533) validated for WB, IHC-P, ELISA in human, mouse, rat. Recombinant protein immunogen, 52kDa. Ideal for keratinocyte and epithelial biology research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 14 (KRT14) Rabbit mAb
Remarks/Alias	K14; NFJ; CK14; EBS1; EBS3; EBS4; EBS1A; EBS1B; EBS1C; EBS1D; Cytokeratin 14 (KRT14)
Species	Human
Gene ID (Human)	3861
Gene ID	3861
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 264-457 of human Cytokeratin 14 (KRT14) (P02533)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02533
Protein Weight	52kDa
Shipping	Ice bag