

Anti-Cytokeratin 16 (Krt16) Monoclonal Antibody For Wb - P08779

Item Number: CM0010857

Introduction

Rabbit monoclonal antibody against human Cytokeratin 16 (KRT16). Validated for WB, IHC-P, IF-P, and ELISA. Cross-reacts with mouse. Suitable for skin biology and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 16 (KRT16) Rabbit mAb
Remarks/Alias	K16; PC1; CK16; K1CP; NEPPK; FNEPPK; KRT16A; Cytokeratin 16 (KRT16)
Species	Human
Gene ID (Human)	3868
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cytokeratin 16 (KRT16) (P08779).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08779
Protein Weight	51kDa
Shipping	Ice bag