

Anti-Cytokeratin 16 (Krt16) Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P - P08779

Item Number: CM0028553

Introduction

Rabbit polyclonal antibody targeting human Cytokeratin 16 (KRT16), suitable for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. For skin biology and keratin research. (SWISS P08779).

[Learn More](#)

Parameter	Details
Product Name	Cytokeratin 16 (KRT16) Rabbit pAb
Remarks/Alias	K16; PC1; CK16; K1CP; NEPPK; FNEPPK; KRT16A; Cytokeratin 16 (KRT16)
Species	Human
GeneID	3868
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-430 of human Cytokeratin 16 (KRT16) (NP_005548.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08779
Protein Weight	51kDa
Shipping	Ice bag