

Anti-Cytokeratin 17 (Krt17) Monoclonal Antibody For Wb, Elisa - Q04695

Item Number: CM0010816

Introduction

Anti-Cytokeratin 17 (KRT17) rabbit monoclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for research on epithelial differentiation, hair follicle biology, and psoriasis.

[Learn More](#)

Parameter	Specification
Product Name	Cytokeratin 17 (KRT17) Rabbit mAb
Remarks/Alias	PC; K17; PC2; 39.1; CK-17; PCHC1; Cytokeratin 17 (KRT17)
Species	Human
GeneID (Human)	3872
GeneID	3872
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cytokeratin 17 (KRT17) (Q04695).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04695
Protein Weight	48kDa
Shipping	Ice bag