

Anti-Cytokeratin 17 (Krt17) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q04695

Item Number: CM0020238

Introduction

Rabbit polyclonal antibody against human Cytokeratin 17 (KRT17), validated for WB, IF/ICC, ELISA. Suitable for studying keratinocyte differentiation, hair follicle biology, psoriasis, and epithelial stem cell research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 17 (KRT17) Rabbit pAb
Remarks/Alias	PC, K17, PC2, 39.1, CK-17, PCHC1, Cytokeratin 17 (KRT17)
Species	Human
GeneID (Human)	3872
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 133-432 of human Cytokeratin 17 (NP_000413.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q04695
Protein Weight	48kDa
Shipping	Ice bag