

Anti-Cytokeratin 17 Monoclonal Antibody For Wb, If, Ihc, Elisa - Q04695

Item Number: CM0008231

Introduction

High-quality rabbit monoclonal antibody against human Cytokeratin 17, validated for WB, IF, IHC, and ELISA. Cross-reacts with mouse and rat. Ideal for epithelial cell research and skin appendage studies.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 17 Rabbit mAb
Remarks/Alias	PC; K17; PC2; 39.1; CK-17; PCHC1
Species	Human
GenelD (Human)	3872
GenelD	3872
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 333-432 of human Cytokeratin 17 (NP_000413.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04695
Protein Weight	48kDa
Shipping	Ice bag