

Anti-Cytokeratin 18 (Ck18) Rabbit Monoclonal Antibody (Wb, If/Icc, Ihc-P, Elisa) - P05783

Item Number: CM0008211

Introduction

Cytokeratin 18 (CK18) rabbit monoclonal antibody validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Targets 48kDa intermediate filament protein in simple epithelia. Suitable for cancer, cell biology, and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 18 Rabbit mAb
Remarks/Alias	K18; CK-18; CYK18
Species	Human
GeneID (Human)	3875
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 241-430 of human Cytokeratin 18 (NP_000215.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05783
Protein Weight	48kDa
Shipping	Ice bag