


Anti-Cytokeratin 19 (Ck19) Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P08727

Item Number: CM0003350



Introduction

Rabbit monoclonal antibody targeting human Cytokeratin-19 (CK19), immunogen synthetic peptide within amino acids 1-100 of KRT19. Validated for WB, IF, IHC-P, ELISA. Suitable for epithelial tissue research.

[Learn More](#)

Parameter	Specification
Product Name	Cytokeratin 19 (CK19) Rabbit mAb
Remarks/Alias	K19; KRT19; K1CS; Cytokeratin 19 (CK19)
Species	Human
GeneID (human)	3880
GeneID	3880
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cytokeratin 19 (KRT19) (KRT19) (P08727).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P08727
Protein Weight	44 kDa
Shipping	Ice bag