


Anti-Cytokeratin 19 (Krt19) Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P08727

Item Number: CM0003427



Introduction

Recombinant rabbit monoclonal antibody targeting human Cytokeratin 19 (KRT19). Suitable for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Recognizes human, mouse, and rat KRT19. Use for intermediate filament and muscle organization research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 19 (KRT19) Rabbit PolymAb®
Remarks/Alias	K19; CK19; K1CS
Species	Human
GeneID (Human)	3880
GeneID	3880
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 243-400 of human Cytokeratin 19 (KRT19) (NP_002267.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08727
Protein Weight	44kDa
Shipping	Ice bag