


Anti-Cytokeratin 19 (Krt19) Rabbit Monoclonal Antibody For If, Ihc, Fc, Elisa - P08727

Item Number: CM0006224



Introduction

Rabbit monoclonal antibody against human Cytokeratin 19 (KRT19), validated for IF/ICC, IHC-P, FC, and ELISA. Reacts with human, mouse, and rat. Ideal for muscle costamere and epithelial biology research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 19 (KRT19) Rabbit mAb
Remarks/Alias	K19; CK19; K1CS
Species	Human
GeneID	3880
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 301-400 of human Cytokeratin 19 (KRT19) (NP_002267.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, IF-P, IHC-P, FC (intra), ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08727
Protein Weight	44kDa
Shipping	Ice bag