

Anti-Cytokeratin 2E (Krt2) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P35908

Item Number: CM0005822



Introduction

Rabbit monoclonal antibody targeting human Cytokeratin 2e (KRT2, CK-2e). Validated in WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for epidermal differentiation and keratinocyte research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 2e (KRT2) Rabbit mAb
Remarks/Alias	K2e; KRTE; CK-2e; KRT2A; KRT2E; Cytokeratin 2e (KRT2)
Species	Human
GenelD	3849
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 293-422 of human KRT2 (P35908)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35908
Protein Weight	65kDa
Shipping	Ice bag