

Anti-Cytokeratin 4 (Krt4) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P19013

Item Number: CM0002007

Introduction

Rabbit monoclonal antibody against human Cytokeratin 4 (KRT4), suitable for WB, IF-P, IHC-P, ELISA. Cross-reacts with human and mouse. Ideal for mucosal stratified squamous epithelium studies.

[Learn More](#)

Parameter	Specification
Product Name	Cytokeratin 4 (KRT4) Rabbit mAb
Remarks/Alias	K4; CK4; CK-4; CYK4; WSN1; Cytokeratin 4 (KRT4)
Species	Human
Gene ID (Human)	3851
Gene ID	3851
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19013
Protein Weight	56kDa
Shipping	Ice bag