


Anti-Cytokeratin 5 (Krt5) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P13647

Item Number: CM0026062



Introduction

Rabbit polyclonal antibody against human Cytokeratin 5 (KRT5), validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse. Ideal for studying intermediate filament biology in basal keratinocytes.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 5 Rabbit pAb
Remarks/Alias	K5; CK5; DDD; DDD1; EBS1; EBS2; EBS2A; EBS2B; EBS2C; EBS2D; EBS2E; EBS2F; KRT5A; Cytokeratin 5 (KRT5)
Species	Human
Gene ID (Human)	3852
Gene ID	3852
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cytokeratin 5 (KRT5) (NP_000415.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13647
Protein Weight	62kDa
Shipping	Ice bag