

Anti-Cytokeratin 6 (Krt6) Monoclonal Antibody For Wb, Ihc-P, Elisa - P48668

Item Number: CM0002536

Introduction

CamelBio's recombinant rabbit monoclonal antibody targets human Cytokeratin 6 (KRT6, SWISS P48668), validated for western blot, immunohistochemistry, and ELISA. Suitable for epithelial biology and dermatology studies.

[Learn More](#)

Parameter	Specification
Product Name	Cytokeratin 6 (KRT6) Rabbit mAb
Remarks/Alias	K6E; KRT6E; PPKNEFD; Cytokeratin 6 (KRT6)
Species	Human
Gene ID (Human)	286887
Gene ID	286887
Immunogen	Synthetic peptide corresponding to amino acids 465-564 of human Cytokeratin 6 (P48668)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P48668
Protein Weight	60 kDa
Shipping	Ice bag