

Anti-Cytokeratin 6 (Krt6) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P02538

Item Number: CM0008141



Introduction

High-purity rabbit monoclonal antibody targeting human Cytokeratin 6 (KRT6A). Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for research on epithelial integrity, wound healing, and keratinocyte biology.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 6 (KRT6) Rabbit mAb
Remarks/Alias	K6A; K6C; K6D; PC3; CK6A; CK6C; CK6D; CK-6C; CK-6E; KRT6C; KRT6D; Cytokeratin 6 (KRT6)
Species	Human
GeneID (Human)	3853
GeneID	3853
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cytokeratin 6 (KRT6) (P02538).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02538
Protein Weight	60kDa
Shipping	Ice bag